

# OCUVANE<sup>®</sup> PLUS

EYELID WIPES

MEDICAL DEVICE

Burning, red, irritated, dry eyes are often caused by inflammation of the eyelids, which can be caused by dirt and allergens deposited there or by a Demodex infection. For lasting alleviation of symptoms, the treatment of eyelid inflammation is essential in order to restore function to the tear glands in the eyelids. Through its special ingredients, **OCUVANE<sup>®</sup> PLUS** eyelid wipes effectively remove dirt and pollen deposited on the edge of the eyelid. Daily use eliminates the blockage of the tear glands in the eyelid, restoring healthy tear production. The active ingredient of the tea tree oil has a disinfectant, anti-inflammatory and regenerative effect that helps to alleviate long-term symptoms.

It can also be used by contact lens wearers to reduce the discomfort associated with wearing lens. With daily treatment, the most common symptoms—such as redness, dryness of the eyes, itching, sensation of tiredness or foreign body—typically decrease within a few days.

## Applications

- To reduce dry, burning sensations in the eyes
- For strained, red eye complaints
- For Demodex infection of the eyelids

## Ingredients

Disodium cocoamphodiacetate, terpinen 4-ol (1%), tocopheryl-acetate (vitamin E), polysorbate 20, citric acid, parfum, purified water

## Instructions for use

Before use, wash your hands, then open the sachet and remove the impregnated **OCUVANE<sup>®</sup> PLUS** eyelid wipes.

### 1. To cleanse eyelids and eyelashes

As shown in the figure, wrap the wipe around your index finger; then, with long smooth strokes, gently wipe the surface of the closed eyelids from the inner to the outer corners.



### 2. Cleansing the edge of the eyelid



In front of a mirror, with fingers placed on the lower eyelash, invert and hold open the lower eyelid. With the wipe wrapped around your index finger, gently wipe the edge of the eyelid in long strokes, several times. Be careful not to touch the eyeball. Using a fresh wipe, also treat the other eye.

In order to remove accumulated dirt from the edges of the eyelids, in mild cases, apply once daily in the evening; in more severe cases, apply twice daily in the morning and evening. It is normal to feel a mild burning sensation during the application of the wipes.

## ⚠ Warnings and precautions

Do not use the wipe after the expiry date that is stated on the packaging, or if the packaging is damaged. Very rarely, hypersensitive reactions to the product's ingredients may occur, with symptoms such as a burning sensation, itching, tearing, redness and swelling of the eyelids. In such cases, use of the wipe should be discontinued. The symptoms outlined above will disappear within a short period of time, but repeated use of the wipe is not recommended in such cases. If, despite treatment, you notice increased inflammatory symptoms (red eyes, the sensation of the presence of a foreign body, tearing, rheumy discharge, pain, or blurred vision) consult an ophthalmologist immediately. Can be used during pregnancy and lactation.

**Store** at room temperature. Keep out of the reach and sight of small children.



## SolMed Pharma Kft.

19 Ady Endre Street, 1024 Budapest  
Hungary  
Email: [info@solmedpharma.hu](mailto:info@solmedpharma.hu)



[www.ocuvane.com](http://www.ocuvane.com)

